



MEVS pro
WIRELESS • SERVICES • INNOVATION

NEW TECHNOLOGY
NO INVESTMENT
UNLIMITED WARRANTY
WIRELESS
MULTIUSER
COMPATIBLE

MANUFACTURED BY:


medicalcables

Innovation and Compatibility for Health



Wiemspro is the only **Whole Body Electrical Muscle Stimulation (EMS)** wireless device, designed and manufactured using state-of-the-art technology and controlled by an Ipad application enabling the trainer to have sessions with groups of up to 12 people at a time, carrying out any type of activity in any location.





The first and only electro-stimulator in the world to incorporate 12 channels, wireless technology and different types of waves (EMS and TENS), all in a single device.

WIEMS pro
WIRELESS • SERVICES • INNOVATION

In the past cables and wires were necessary to connect the user to the machine and this hindered freedom of movement. Now Wiemspro has eliminated this by means of the use of Bluetooth technology.

wireless



The device has an exclusive design using a small latest generation lithium-ion battery that offers 7 hours of continued use. This lets the trainer carry out sessions indoors in a gym or in the home or outdoors on the beach, in the countryside, etc.

An application has been developed for Ipad which enables the device to be controlled in a simple and intuitive way. Thanks to this application and the device, a trainer can work with up to 12 users all at the same time, enabling group and individual sessions to be directed from indoors or outdoors.



outdoor

multi-user



Compatibility

Wiemspro comes with a series of cables (accessories) that make it compatible with the suits manufactured by major brands already in the market. This enables the trainer to use suits already owned regardless of the brand.



Ipad Control

Control is done by means of a powerful and easy to use multi-device control interface. The Ipad application offers, at a single glance, complete information of the basic training parameters of each person being trained and this enables quick and simple readjustments.



EMS / TENS

Wiemspro makes it possible to have EMS and TENS in the same device, offering a high quality wave and unprecedented sensations with this type of device, due to the utilization of the latest technology in signal generation and processing.

programs

Wiemspro has 20 training programs, and additional programs that are fully configurable by the trainer.

WIEMS pro
WIRELESS • SERVICES • INNOVATION





micro-usb

Universal charging system with micro-USB connector.

easy

Compact, attractive and easy to use.

certifications

Complies with all required certifications regulated by the EC (Europe).



Wiemspro comes with a series of innovations that make it a unique device in the sector for body electro-stimulation and it allows the trainer to carry out single and group sessions all at the same time, either indoors or outdoors, enabling the trainer to use already owned vests regardless of the brand.



CONTENT

- 1 wiemspro device.
- 1 AC adaptor with micro-USB connector.
- User guide.

ACCESSORIES

Compatible cables
for different suits.

TECHNICAL INFORMATION

- Wireless data transfer using Bluetooth
- 12 channels.
- EMS + TENS technology in a single device.

WARRANTY

Unlimited for the duration of the contract

WIEMS pro

WIRELESS • SERVICES • INNOVATION



innovation@wiemspro.com · www.wiemspro.com

medicalcables
Innovation and Compatibility for Health

